

## Beer's Law calculations

This is a review of the steps used in solution calculations and using Beer's Law to analyze a  $\text{K}_2\text{CrO}_4$  solution of unknown molarity. First one prepares a series of standard  $\text{K}_2\text{CrO}_4$  solutions.

1. Preparation of 50.00 mL of 0.400 M  $\text{K}_2\text{CrO}_4$  solution.

$$(0.400\text{mol salt/L})(0.0500\text{ L})(194\text{g salt/mol}) = 3.88\text{ g K}_2\text{CrO}_4\text{ needed.}$$

Add 3.88 g of the salt to a 50.00 mL volumetric flask, add enough water to dissolve it and then dilute to the calibration mark with water. Mix well.

2. Prepare 50.00 mL of 0.300 M solution of the salt from the 0.400 M solution.

Use the dilution formula:  $M_1V_1 = M_2V_2$

$$(0.400\text{M})(V_1) = (0.300\text{M})(50.00\text{mL})$$

$V_1 = 37.5\text{ mL}$  (the volume of the 0.400M salt which must be diluted to 50.00 mL to give a 0.300M solution. 37.5 mL of the 0.400M solution is added to a 50.00 mL volumetric flask and water added to the calibration mark.

3. Preparation of a calibration curve.

Using the techniques in steps 1 and 2 above a series of different concentrations of the salt are made and their absorbance measured to give the following data.

$[\text{K}_2\text{CrO}_4], \text{M}$	Absorbance, A
0.000	0.000
0.100	0.145
0.200	0.255
0.300	0.415
0.400	0.525

A graph of Absorbance (y-axis) versus concentration (x-axis) is then plotted and its slope,  $\Delta y/\Delta x$ , is evaluated. See **Using Graphical Analysis 3.1.1** or **Graphing in Excel** references for plotting the graph. Graphs using these two programs are shown on pages 2 and 3.


#### 4. Evaluation of the molarity of an unknown $K_2CrO_4$ solution.

The absorbance of the unknown solution is measured under the same conditions as the standard solutions. Its absorbance is 0.250. Since slope ( $m$ ) = Absorbance / concentration,  $[K_2CrO_4] = \text{absorbance}/\text{slope} = 0.250/1.32/M = 0.189M$

A more accurate method is using the  $y = mx + b$  formula obtained from the plotted graph where  $y$  is absorbance and  $x$  is the concentration.

Thus  $0.250 = 1.32x + 0.004$  and  $x$  is 0.186 M.

### Beer's Law Using Graphical Analysis 3.1.1



## Beer's Law Using Excel

